



# OMNISTAR™/ THERMOSTAR™

The efficient solution for gas analysis.  
Intelligent software. Wide range of applications.



**The efficient solution for gas analysis.**  
**Intelligent software. Wide range of applications.**

## **What are OmniStar™ and ThermoStar™?**

OmniStar and ThermoStar are compact benchtop analysis systems for sample gases at, or below atmospheric pressure. They are the perfect complete solution for gas analysis, especially in chemical processes, the semiconductor industry, in metallurgy, fermentation, catalysis, laser technology and environmental analysis.

The analysis systems consist of an inlet system, a PrismaPlus mass spectrometer, a dry-compressing diaphragm vacuum pump and a HiPace turbopump.

The gas inlet is equipped with a capillary that can be heated to 350°C; this capillary is made of stainless steel in the OmniStar and of quartz glass in the ThermoStar. The heated capillary prevents vapors from condensing during analysis of the sample gas.

This 2-stage inlet enables a quasi segregation-free gas supply. A pressure-regulated gas inlet is optionally available for varying inlet pressures from 1000 hPa (mbar) down to  $10^{-3}$  hPa (mbar).

The Quadera mass spectrometer software enables both qualitative and quantitative analyses to be performed on up to 128 different masses. The systems cover the mass ranges of 1-100 amu, 1-200 amu and 1-300 amu

## **ThermoStar™ – Special solution**

The ThermoStar solution was specifically developed for coupling with thermo balances. The inlet system incorporates a quartz capillary and platinum orifice, assuring that even minute concentrations can still be analyzed. In contrast to competitive analysis methods, such as FTIR and IR, it is possible to detect all gases within the mass range.



Laser technology



Catalysis



Semiconductor industry



#### Advantages at a glance

- Qualitative and quantitative gas analysis, connection port for calibration gas inlet systems
- Low detection limit ( $< 1$  ppm), even for condensable gases
- Compact, easy-to-operate analysis unit
- Heated capillary inlet, up to  $350^{\circ}\text{C}$
- Bakeable all-metal sealed high vacuum chamber for low backgrounds
- Utmost sensitivity thanks to enclosed ion source and field axis technology
- Reliable identification of unknown gases
- Monitoring of up to 128 masses
- Mass ranges of 1 to 100 amu, 1 to 200 amu and 1 to 300 amu
- Applications from vacuum up to atmospheric pressure
- Gas analysis under varying sample pressures – pressure-regulated gas inlet (optional)
- Fast, reliable and precise measurement of non-polar and noble gases
- Soft Ionization option


## Analysis software for total analysis

### Quadera® analysis software

The Quadera analysis software enables you to easily monitor your processes on a PC in realtime. The unit's short response time enables the sample analysis to be quickly conducted. With its modular design, the Quadera software offers a clearly organized, user-friendly platform for the capture and presentation of measurement data and parameter records. Complete measurement procedures can be programmed.




#### Xe in air/OmniStar™



The graph (mass range 125–140 amu) shows the xenon isotopes and the excellent detection limit of the OmniStar. Xenon has stable isotopes at 129, 131, 132, 134 and 136 amu. The concentration of <sup>136</sup>Xe in air is 7.8 ppb.

#### Respiratory Gas Analysis/OmniStar™



Three breathing cycles of an endurance test. The breath is characterized by the gases O<sub>2</sub> and CO<sub>2</sub>. The argon in the air serves as a reference gas for demonstrating the stability of the system. Very quickly measured by the OmniStar at atmospheric pressure in the Multiple Concentration Determination (MCD) mode.

### Advantages at a glance

- User-friendly, intuitive operation
- Customer-specific user interface
- Automated measurement procedures, programmable by means of the integral Microsoft® VSTA® script editor
- I/O module for exchanging measurement data
- Integral library of mass spectral data
- Simple definition of measurement recipes
- Mass spectrometer data can be linked with external signals

# THERMOSTAR™


The special solution in connection with thermo balances.

## ThermoStar™

The ThermoStar gas analysis unit is a version that is specially designed for being coupled with thermo balances. Hightemperature gas samples can be admitted by means of a quartz capillary.



## TGA-MS




The gases water vapor (18), carbon monoxide (28) and carbon dioxide (44) that evolve when conducting a thermal analysis of calcium oxalate are shown as a function of time together with weight and temperature.

## Advantages at a glance

- Even small concentrations of reactive and condensable gases can be detected
- Inert inlet, no change in gas composition
- Gas inlet heatable up to 350 °C
- Multi-gas analysis

The efficient solution for gas analysis.  
Intelligent software. Wide range of applications.

## Dimensions



Dimensions in mm

## Technical data

Gas analysis system	OmniStar™	ThermoStar™
Mass ranges, amu	1 – 100 / 1 – 200 / 1 – 300	
Gas connection	Stainless steel capillary	Quartz capillary
Gas inlet	Via software-controlled inlet valve, or user interface	Continuously open
Pressure reduction	2-stage, segregation-free	
Gas flow rate, sccm	1 – 2	
Sample gas pressure, hPa (mbar)	Up to 1,000	
Capillary operating temperature, °C	Up to 350	
Analyzer	QMA 200 M gas tight	
Rod system, material / diameter / length, mm	Stainless steel/6/100	
Detector	C-SEM/Faraday	
Mass spectrometer electronics	QME 220 M	
Software	Quadera®	
Contribution to neighboring mass: 40 to 41	<10 ppm / <20 ppm / <50 ppm	
Min. detection limit, C-SEM	<1 ppm / <1 ppm / <1 ppm	
Min. detection limit, Faraday	<20 ppm / <40 ppm / <100 ppm	
Resolution, adjustable at 10 % peak height, amu	0.5 – 2.5	
Dimensions (L x B x H), mm	545 x 305 x 395	
Weight, kg	35	
Mains requirement: voltage (range), VAC	100 – 230	
Interface	Ethernet	
	Analog input: 5x ± 10 V / 12 bit	
	Analog output: 4x 0...10 V / 12 bit	
	Digital input: 4x	
	Digital output: 7x sink, optical insulated, 24 V	

System	a
OmniStar	8
ThermoStar	9
Analyzer/Filament	bc
Filament: Yttriated iridium, integrated option for mass calibration	01
Filament: Yttriated iridium	02
Filament: Tungsten, integrated option for mass calibration	05
Filament: Tungsten	06
Filament: Yttriated iridium, integrated option for mass calibration, corrosive gas version with inert gas inlet	11
Filament: Yttriated iridium, corrosive gas version with inert gas inlet	12
Filament: Tungsten, integrated option for mass calibration, corrosive gas version with inert gas inlet	15
Filament: Tungsten, corrosive gas version with inert gas inlet	16
Gas inlet system	def
ThermoStar, quartz, 1 m temperature-regulated gas sampling line, 200 °C	111
ThermoStar, quartz, 2 m temperature-regulated gas sampling line, 200 °C	112
ThermoStar, quartz, 1 m temperature-regulated gas sampling line, 350 °C	113
OmniStar, quartz (special), 1 m temperature-regulated gas sampling line, 200 °C	151
OmniStar, quartz (special), 2 m temperature-regulated gas sampling line, 200 °C	152
OmniStar, stainless steel, 1 m temperature-regulated gas sampling line, 200 °C	171
OmniStar, stainless steel, 2 m temperature-regulated gas sampling line, 200 °C	172
OmniStar, stainless steel, 1 m temperature-regulated gas sampling line, 350 °C	173
Pressure-regulated gas inlet, 5 decades, $5 \cdot 10^{-3}$ mbar up to atmospheric pressure	555
Mass range	g
1–100 amu	1
1–200 amu	2
1–300 amu	3
Front panel	h
OmniStar	1
ThermoStar	2



OmniStar™ and ThermoStar™ are trademarks of INFICON  
 Quadera® – An INFICON product  
 VSTA® – A Microsoft® product





伯东企业（上海）有限公司  
HAKUTO ENTERPRISES (SHANGHAI) LTD.  
TEL: (86) 21-5046 3511 FAX: (86) 21-5046 1490

## 上海伯东—您理想的真空合作伙伴！

伯东公司成立于1953年，全球范围内有一千七百名员工。伯东企业（上海）有限公司 [www.hakuto-vacuum.cn](http://www.hakuto-vacuum.cn) 主要从事半导体设备、电子零部件加工、销售及相关服务工作，集团年营业额约为100亿人民币。伯东公司在国内主要城市：上海、北京、苏州、广州、深圳、成都、厦门及香港、台湾设有营业据点。

上海伯东真空产品事业部秉承为广大中国客户提供世界一流的真空产品、推动中国真空工艺发展、承担企业社会责任为经营理念，已累计为超过10,000家企业提供真空服务，覆盖工业、半导体、镀膜、科研和分析行业。

上海伯东 Hakuto 代理品牌包含：德国 Pfeiffer 全系列真空产品、美国 KRI 离子源、美国 Brooks Polycold 制冷机、美国 HVA 真空阀门、美国 inTEST (Temptronic) 高低温循环测试机，NS 离子蚀刻机，真空镀膜系统及配套真空零配件等。

上海伯东 Hakuto 作为国际知名真空品牌指定授权代理商，全权负责在中国地区的销售和维修服务，100%使用原厂进口零部件以及国外受训维修工程师；拥有完全的拆解、维修、校正能力并提供24小时在线服务。







伯东企业（上海）有限公司  
**HAKUTO ENTERPRISES (SHANGHAI) LTD.**  
TEL: (86) 21-5046 3511      FAX: (86) 21-5046 1490

**伯东企业（上海）有限公司 销售及售后服务中心**

上海市外高桥保税区希雅路33号17#3楼B座

电 话: +86-21-5046-3511

传 真: +86-21-5046-1490

邮箱: [ec@hakuto-vacuum.cn](mailto:ec@hakuto-vacuum.cn)

官 网: <http://www.hakuto-vacuum.cn>

24 h售后维修专线: 1381-880-1372

**伯东企业（上海）有限公司 北京办事处 销售中心**

北京丰台区芳群园四区22号楼金城中心705室

电 话: +86- 010-8768-4988

传 真: +86-21-5046-1490

邮箱: [ec@hakuto-vacuum.cn](mailto:ec@hakuto-vacuum.cn)

官 网: <http://www.hakuto-vacuum.cn>

**伯东企业（上海）有限公司 深圳办事处 销售中心**

深圳市罗湖区宝安南路1036号鼎丰大厦501

电 话: +86- 0755-2547-3988. 643

传 真: +86-21-5046-1490

邮箱: [ec@hakuto-vacuum.cn](mailto:ec@hakuto-vacuum.cn)

官 网: <http://www.hakuto-vacuum.cn>

**伯东国际通商股份有限公司 销售及售后服务中心**

台北市中山区八德路二段308号6楼

电 话: +886-2-8772-8910

传 真: +886-2-8772-8918

邮箱: [ec@hakuto-vacuum.cn](mailto:ec@hakuto-vacuum.cn)

官 网: <http://www.hakuto-vacuum.com.tw>

新竹市埔顶路25号8楼之一

电 话: +886-3-5679508

传 真: +886-3-5670049

邮箱: [ec@hakuto-vacuum.cn](mailto:ec@hakuto-vacuum.cn)

官 网: <http://www.hakuto-vacuum.com.tw>